

NE-300 "Just Infusion"TM Syringe Pump



NE-300 "Just Infusion"TM Syringe Pump Pump NE-300

International orders please contact us for more information

Description:

"Just the features that you want, and none that you don't want to pay for."

- Holds small or big syringes up to 60 mL
- Pump continuously until you stop the pump
- Infusion rate can be changed while pumping
- Remembers previous settings on power up
- Power failure mode allows pump to continue pumping after a reset
- Syringe purge mode
- Displays total volume dispensed in mL or μ L units
- Maximum pumping rate of 1500 mL/hr with a 60 mL syringe
- Precise reproducible flow rates
- Selectable infusion rates units: mL/hr, μ L/hr, mL/min, μ L/min
- Pictured Syringe is for Demonstration purposes only
- **"Infuses Only!"**
- -Does not withdraw, volume target, or have computer interface
- -Just the features that you want to pay for
- -See Syringe and Plumbing supplies to purchase syringes.

Fabrizio Barbieri (ceo)

KF Technology Srl

T.(+39) 06 454.34.179 - info mobile: (+39) 339 533.03.22 -

fax.(+39) 06 9725.3131- SkyPe : kft2002

fabrizio@kftechnology.it - www.kftechnology.it