

Rat Standard Chambers



Front / Top Loading Operant Chamber Model 80003

Built for long-lasting durability and easy cleaning, our front / top-loading operant chamber is constructed of .032 in stainless steel and .25 in polycarbonate. The base is made of high density polyethylene that is stain- and scent-resistant. Comes included with removable shockable grid rod floor.



Front / Top Loading Operant Chamber for Rats with No Shock Floor Model 80003NS

Built for long-lasting durability and easy cleaning, our front / top-loading operant chamber is constructed of .032 in stainless steel and .25 in polycarbonate. The base is made of high density polyethylene that is stain- and scent-resistant. Comes included with No Shock Floor.

Fabrizio Barbieri (ceo)

KF Technology Srl

Office: Via Amedeo Bocchi 84 - piano terra - int. 8 - 00125 Roma

(+39) 06 454.34.179 - info mobile:(+39) 339 533.03.22 - fax.(+39) 069725.3131- SkyPe : kft2002

fabrizio@kftechnology.it - www.kftechnology.it

IMPORTANTE: è disponibile il ns catalogo on-line..... www.kftechnology.it